Lucitone® HIPA

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Mode d'emploi

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INDICATIONS:

Lucitone HIPA Denture Base is a self-curing denture base material designed for use in fabrication, repair, rebasing, or relining of full and partial dentures including implant overdentures or other dental appliances.

CONTRAINDICATIONS:

HIPA Denture Base is contraindicated for patients and users who have a history of allergic reaction to methyl methacrylate monomer.

WARNINGS:

1. Material contains polymerizable monomers which may cause skin sensitization (allergic contact dermatitis) or other allergic reactions in susceptible persons. Wear protective gloves and protective clothing. Wash thoroughly with soap and water after contact. If dermatitis or other symptoms persist, seek medical advice.
2. Avoid inhalation or ingestion. High vapor concentration can induce headache, nausea, and irritation of eyes and respiratory system. Work in a well ventilated area. Liquid contact with eyes may cause possible corneal damage. Wear eye protection. Excessive long-term exposure may be associated with other more serious health effects.

Inhalation: Move subject to fresh air. Give oxygen or artificial respiration as required.

Ingestion: Contact your regional Poison Control Center immediately.

Eye Contact: Flush eyes promptly with copious amounts of water for 15 minutes, and consult a physician. Wash skin with soap and water.

PRECAUTIONS:

1. When grinding prosthodontic resins, proper ventilation, masks, and vacuum systems should be used.
2. Store at 60°F-80°F (16°C-27°C), away from moisture and direct sunlight. The liquid contains unsaturated monomers which may polymerize prematurely if stored at excessively high temperature and/or in sunlight. The powder contains an organic polymerization initiator which may degrade if stored at excessively high temperature.
3. Warm materials (powder or liquid) will reduce working/pouring time significantly. When pouring, the best results are achieved when material temperatures are at or below 73°F/ 23°C.
4. For disposal, please follow the instructions on the safety data sheet or national regulations.

- Liquid contains methyl methacrylate monomer, a flammable liquid with a flash point of 50°F (10°C). Keep away from heat, sparks, and open flame.
- Use liquid in a well-ventilated area.
- Replace liquid bottle cap when not in use.
- Excess liquid already poured out should not be returned to the original bottle or slight discoloration may occur.
- Use entire contents of bottle within 5 months of opening or slight discoloration may occur.
- To ensure uniformity shake powder jar (with a rocking motion) to prevent concentration of small particles at the bottom of the jar.
- DO NOT** use detergents to clear residual wax from casts. Detergent residue may impact the color uniformity.
- Allergic contact dermatitis and other allergic reactions may occur in susceptible individuals. Residual monomer in fully cured materials can be minimized by soaking the cured prosthesis in warm water for several days.
- Particulates will be generated when grinding acrylate resins. Eye, skin and respiratory irritation may occur if appropriate engineering controls are not used.

STEP BY STEP INSTRUCTIONS:

PROPERTIES AND RECOMMENDED PROCEDURES

The following table describes typical properties of Lucitone HIPA Denture Base, plus recommended temperatures and other important processing information. Further detail of this information is found later in the text.

Maximum Residual MMA Content	< 4.5% (% mass fraction)
Type and Class (according to ISO 20795-1)	Type 2, Class 2
Storage temperature for powder and liquid	60°F-80°F (16°C-27°C)
Powder / Liquid ratio (HIPA Denture Base Powder and Liquid)	20 g (28 cc) / 15 ml

ADVERSE REACTIONS:

- Corneal damage, headache, nausea, vomiting may occur with exposure to methyl methacrylate monomer. (See Warnings and Precautions statements.)

Mixing time (time required to thoroughly mix)	15 seconds
Time to reach pourable stage @ 73°F ± 2°F (23 ± 1°C)	can be poured immediately
Pour Time	about 3 minutes @ 73°F ± 2°F (23 ± 1°C)
Material used to prepare mold	hydrocolloid or silicone
Temperature of hydrocolloid when pouring	60°F-80°F (16°C-27°C)
Cure time, temperature, and pressure (HIPA Denture Base)	pressure pot, 30 minutes Minimum of 20 psi (1.4 bars) 113°F ± 2°F (45°C ± 1°C) or 131°F ± 2°F (55°C ± 1°C) if using pour-type duplicating silicone

FABRICATION OF FULL, PARTIAL AND IMPLANT OVER DENTURES - POURING PROCESS

HYDROCOLLOID / SILICONE INVESTMENT AND SPRUING

1. Incorporate a slight taper into the base of the cast to facilitate removal from the mold material. Taper so the base is the widest.
2. Soak the waxed stone cast in water for 10-20 minutes prior to pouring the mold material (hydrocolloid or silicone).
3. Follow the manufacturer's instructions for preparation and use of the selected hydrocolloid or silicone product and flask.
 - a. Create or place sprue holes in the hydrocolloid or silicone mold as shown in Figure 1.
 - b. Sprue diameter should be large to reduce tendency for air to be entrapped, for example 7-8 mm where possible.

IMPORTANT: When fabricating partials each acrylic area must have two sprues to avoid trapped air.


Figure 1. Spruing Techniques

WAX ELIMINATION AND SEPARATOR APPLICATION

1. Boil-off wax with clean boiling water. DO NOT use detergent to clean casts. Detergent residue may impact color uniformity and could cause bleaching on tissue side of the denture.
2. Clean teeth with clean boiling water. Care must be taken to assure that all wax is removed from the teeth.
3. Add mechanical retention to the teeth. Remove only enough tooth material to enhance retention. Excessive removal of tooth material will weaken the tooth and may lead to fracture.
4. Immediately after boil out, place the cast in a bowl of water for 10 minutes to eliminate air; extend soaking time if needed until bubble free.

Note: If there is any air inside the cast, small bubbles may occur in the denture.

5. Remove the cast from the water and allow to dry for 1-2 minutes.
 6. Use a brush to apply a thin coat of AL-COTE[®] Separating Agent to the cast. Do not allow the separator to puddle. Apply a second layer as needed. Make sure the separator is completely dry before placing it back in the mold.
- Note:** Blanching may occur if separator is not completely dry.
7. Blow out the hydrocolloid mold with air until it is completely dry.
 8. Place prepared teeth into mold and place cast back into the hydrocolloid (or silicone) mold making sure the cast is seated properly.
- Note:** Make sure the mold is not warm when pouring the acrylic.
9. If desired, place a plastic sheet over base of model to ensure that acrylic will not leak when being poured.
 10. Re-assemble the flask being especially careful to avoid possible tooth and / or cast movement.

MIXING POWDER AND LIQUID

1. For the fabrication of one full denture, use one (1) Large Scoop (or 20g) of HIPA Denture Base Powder and 15 ml of HIPA Denture Base Liquid measured with the HIPA liquid vial.
2. For the fabrication of a partial denture, use one (1) Small Scoop (or 11g) of HIPA Denture Base Powder and 8 ml of HIPA Denture Base Liquid measured with the HIPA liquid vial.
Note: After measuring the liquid, invert the vial on a paper towel to allow residual liquid to empty from the vial. Do not return excess liquid to original bottle.
3. Add powder to liquid and mix thoroughly for a minimum of 15 seconds to completely wet out the powder. Avoid air bubble entrapment.
Note: Failure to mix well will cause white areas or spots in the denture.

POURING

Pour acrylic into mold (pour into one side sprue hole only). Pour time is about 3 minutes at 23°C (73°F).

CURING

1. Check the water temperature in the pressure pot with an accurate thermometer and adjust to 113°F \pm 2°F (45°C \pm 1°C) or 131°F \pm 2°F (55°C \pm 1°C) when using pour duplicating

silicone.

2. Place the flask in the water with sprue holes upright at 113°F \pm 2°F (45°C \pm 1°C). The water level should be below sprue height but above 2/3 the height of the flask. DO NOT FULLY SUBMERGE THE FLASK IN THE WATER.
3. Immediately close the pressure pot and apply 20 psi (1.4 bars) or greater pressure.
Note: Delay in applying pressure may introduce porosity.
4. Cure for 30 minutes.
5. Cool flask and remove denture from the mold.

FINISHING AND POLISHING

Finish and polish the denture using normal procedures.

REPAIR PROCEDURE

REPAIRING A FRACTURE

PREPARATION OF THE REPAIR AREA

1. To repair a broken denture, hold the fractured parts in position with splints and sticky wax while pouring the stone cast.

2. Apply AL-COTE® Separating Agent to the cast.
3. Open the repair area to approximately 1/4 inch, and grind a bevel on the edges of opening with narrower width on the tissue side of the repair.
4. The contact or surfaces to be repaired must be well roughened.

PRE-MIX DOUGH PROCEDURE

1. Measure one (1) Small Scoop (or 11g) of HIPA Denture Base Powder and 8 ml of HIPA Denture Base Liquid measured with the HIPA liquid vial.
Note: After measuring the liquid, invert the vial on a paper towel to allow residual liquid to empty from the vial.
2. Add powder to liquid and mix thoroughly for a minimum of 15 seconds to completely wet out the powder. Avoid air bubble entrapment.
Note: Failure to mix well will cause white areas or spots in the denture.
3. Allow the P/L mixture to bench set to reach the desired consistency for repairing.
4. Wet the surface of the prepared area with HIPA Denture Base Liquid then place the mixture into the repair area. Build to a slight over contour.

DUSTING (SALT & PEPPER) PROCEDURE

Build up repair in successive layers by wetting the roughened surface with the HIPA Denture Base Liquid and applying a thin layer of HIPA Denture Base Powder. Repeat this layering until the desired thickness is achieved. The repair area should be slightly over contoured for best results.

CURING

1. Check the water temperature in the pressure pot with an accurate thermometer and adjust to 113°F ± 2°F (45°C ± 1°C).
2. Submerge the repaired denture, mounted upon its stone cast, in a pressure pot with water temperature at 113°F ± 2°F (45°C ± 1°C).
3. Immediately close the pressure pot and apply 20 psi (1.4 bars) or greater pressure.
Note: Delay in applying pressure may introduce porosity.
4. Cure for 30 minutes.

FINISHING AND POLISHING

1. Grind off excess of the over-contoured material down to the level of denture surface.
2. Finish and polish the denture using

normal procedures.

REPLACING A TOOTH

PREPARATION OF REPAIR AREA

1. Grind the denture and make a box preparation on the palatal or lingual side of the denture.
2. Select a replacement tooth and add mechanical retention to the tooth. Remove only enough tooth material to enhance retention. Excessive removal of tooth material will weaken the tooth and may lead to fracture.
3. Use sticky wax to position the tooth in place from the lingual.
4. Prepare a plaster cast facial matrix to hold the tooth in position during repair. The matrix should include the tooth that is being repaired and the tooth on each side.
5. Clean all wax from the box preparation and tooth.
6. After the plaster has set, remove the matrix and apply AL-COTE® Separating Agent to the matrix.
7. Reassemble the denture, the tooth, and the matrix in correct alignment. Use sticky wax to hold parts in position.

PRE-MIX DOUGH PROCEDURE

1. Measure one (1) Small Scoop (or 11g) of HIPA Denture Base Powder and 8 ml of HIPA Denture Base Liquid measured with the HIPA liquid vial. **Note:** After measuring the liquid, invert the vial on a paper towel to allow residual liquid to empty from the vial.
2. Add powder to liquid and mix thoroughly for a minimum of 15 seconds to completely wet out the powder. Avoid air bubble entrapment. **Note:** Failure to mix well will cause white areas or spots in the repair.
3. Allow the P/L mixture to bench set to reach the desired consistency for repairing.
4. Wet the surface of the prepared area with HIPA Denture Base Liquid then place the mixture into the repair area. Build to a slight over contour.

DUSTING (SALT & PEPPER) PROCEDURE

Build up repair in successive layers by wetting the roughened surface with the HIPA Denture Base Liquid and applying a thin layer of HIPA Denture Base Powder. Repeat this layering until the desired thickness is achieved. The repair area should be slightly over contoured for best results.

CURING

1. Check the water temperature in the pressure pot with an accurate thermometer and adjust to $113^{\circ}\text{F} \pm 2^{\circ}\text{F}$ ($45^{\circ}\text{C} \pm 1^{\circ}\text{C}$).
2. Submerge the denture with tooth repair in a pressure pot with water temperature at $113^{\circ}\text{F} \pm 2^{\circ}\text{F}$ ($45^{\circ}\text{C} \pm 1^{\circ}\text{C}$).
3. Immediately close the pressure pot and apply 20 psi (1.4 bars) or greater pressure.
Note: Delay in applying pressure may introduce porosity.
4. Cure for 30 minutes.

FINISHING AND POLISHING

1. Grind off excess of the over-contoured material down to the level of denture surface.
2. Finish and polish the denture using normal procedures.

INDIRECT RELINE PROCEDURE

HIPA Denture Base is recommended for relining poorly fitting dentures by the indirect method. A reline jig is recommended.

An impression is obtained, using old denture as tray. Periphery is trimmed and large undercuts eliminated before taking impression.

PREPARATION OF AREA TO BE RELINED

1. Grind the tissue surface to relieve the denture.
2. Take an impression in the relieved denture.
3. Pour a cast.
4. Then mount the cast, with the denture in position, on an articulator or duplicator to preserve the tooth/cast dimension.
5. Make an occlusal index by placing teeth in plaster.
6. Remove the denture from the cast and thoroughly clean out the impression material.
7. All remaining undercuts are then removed.
8. Replace the denture in the prepared occlusal index, close the articulator and seal the denture to the cast. Wax the borders of the denture.

INVESTING AND SEPARATOR APPLICATION

1. Invest the case in the customary manner. Heat the flask in hot water, separate and remove all remaining impression material and wax.

2. Apply AL-COTE Separator to the exposed investment and on the model. Allow the AL-COTE to dry. Bench cool the flask to room temperature.

RELINE

1. Wet the tissue side of the denture with the HIPA Denture Base Liquid.
2. Measure one (1) Small Scoop (or 11g) of HIPA Denture Base Powder and 8 ml of HIPA Denture Base Liquid measured with the HIPA liquid vial.
Note: After measuring the liquid, invert the vial on a paper towel to allow residual liquid to empty from the vial.
3. Add powder to liquid and mix thoroughly for a minimum of 15 seconds to completely wet out the powder. Avoid air bubble entrapment.
Note: Failure to mix well will cause white areas or spots in the reline.
4. Allow P/L mixture to bench set until the desired consistency for relining. When P/L mixture achieves the desired consistency for relining (softer texture than that of an ordinary packable denture resin), place the HIPA acrylic into the denture, reseal the cast / denture on the jig and close the jig.

CURE

1. Check the water temperature in

the pressure pot with an accurate thermometer and adjust to $113^{\circ}\text{F} \pm 2^{\circ}\text{F}$ ($45^{\circ}\text{C} \pm 1^{\circ}\text{C}$).

2. Submerge the relined denture / jig in a pressure pot with water temperature at $113^{\circ}\text{F} \pm 2^{\circ}\text{F}$ ($45^{\circ}\text{C} \pm 1^{\circ}\text{C}$).
3. Immediately close the pressure pot and apply 20 psi (1.4 bars) or greater pressure.
Note: Delay in applying pressure may introduce porosity.
4. Cure for 30 minutes.
5. Cool flask and remove relined denture/ jig.

FINISHING AND POLISHING

Finish and polish the relined denture in usual manner.

REBASE PROCEDURE

HIPA Denture Base is recommended for rebasing poorly fitting dentures by the indirect method.

PREPARATION OF DENTURE TO BE REBASED

1. Make a rebase impression.
2. Prepare a master cast in the same manner as for a complete denture (beading and boxing). Do not remove

- the denture or the impression from the master cast after it is poured.
- Prepare an occlusal stint before the denture rebase is invested in the flask.
 - Remove the master cast and denture from the occlusal stint, then remove the denture from the master cast. Thoroughly clean out the impression material.
 - The old denture base is ground away leaving only the teeth and enough denture base to retain the teeth in one piece.
 - Place the denture teeth into the occlusal stint and secure it with sticky wax.
 - Place master cast and occlusal stint on the articulator or duplicator. Lute the teeth to the master cast.
 - Remove the master cast with teeth luted in position from the stint.
 - Contour the new denture base with pink wax in the normal manner.
 - Flask the waxed "rebased" denture base per directions below.
- Soak the waxed stone cast in water for 10-20 minutes prior to pouring the mold material (hydrocolloid or silicone).
 - Follow the manufacturer's instructions for preparation and use of the selected hydrocolloid or silicone product and flask.
 - Create or place sprue holes in the hydrocolloid or silicone mold as shown in Figure 1.
 - Sprue diameter should be large to reduce tendency for air to be entrapped, for example 7-8 mm where possible.

IMPORTANT: When fabricating partials each acrylic area must have two sprues to avoid trapped air.


Figure 1. Spruing Techniques

HYDROCOLLOID / SILICONE INVESTMENT AND SPRUING

- Incorporate a slight taper into the base of the cast to facilitate removal from the mold material. Taper so the base is the widest.

WAX ELIMINATION AND SEPARATOR APPLICATION

1. Boil-off wax with clean boiling water. DO NOT use detergent to clean casts. Detergent residue may impact color uniformity and could cause bleaching on tissue side of the denture.
2. Immediately after boil out, place the cast in a bowl of water for 10 minutes to eliminate air; extend soaking time if needed until bubble free.
Note: If there is any air inside the cast, small bubbles may occur in the denture.
3. Remove the cast from the water and allow to dry for 1-2 minutes.
4. Use a brush to apply a thin coat of AL-COTE® Separating Agent to the cast. Do not allow the separator to puddle. Apply a second layer as needed. Make sure the separator is completely dry before placing it back in the mold.
Note: Blanching may occur if separator is not completely dry.
5. Blow out the hydrocolloid mold with air until it is completely dry.
6. Place prepared denture teeth/ remaining denture base into the mold and place cast back into the hydrocolloid (or silicone) mold making

sure the cast is seated properly.

Note: Make sure the mold is not warm when pouring the acrylic.

7. If desired place a plastic sheet over base of model to ensure that acrylic will not leak when being poured.
8. Re-assemble the flask being especially careful to avoid possible teeth/denture base and / or cast movement.

MIXING POWDER AND LIQUID

1. For the fabrication of one full denture, use one (1) Large Scoop (or 20g) of HIPA Denture Base Powder and 15 ml of HIPA Denture Base Liquid measured with the HIPA liquid vial.
2. For the fabrication of a partial denture, use one (1) Small Scoop (or 11g) of HIPA Denture Base Powder and 8 ml of HIPA Denture Base Liquid measured with the HIPA liquid vial.
Note: After measuring the liquid, invert the vial on a paper towel to allow residual liquid to empty from the vial. Do not return excess liquid to original bottle.
3. Add powder to liquid and mix thoroughly for a minimum of 15 seconds to completely wet out the powder. Avoid air bubble entrapment.
Note: Failure to mix well will cause

white areas or spots in the denture.

4. Wet the surface of the prepared teeth/denture base with HIPA Denture Base Liquid.

POURING

Pour acrylic into mold (pour into one side sprue hole only). Pour time is about 3 minutes at 23°C (73°F).

CURING

1. Check the water temperature in the pressure pot with an accurate thermometer and adjust to 113°F ± 2°F (45°C ± 1°C) or 131°F ± 2°F (55°C ± 1°C) when using pour duplicating silicone.
2. Place the flask in the water with sprue holes upright at 113°F ± 2°F (45°C ± 1°C). The water level should be below sprue height but above 2/3 the height of the flask. **DO NOT FULLY SUBMERGE THE FLASK IN THE WATER.**
3. Immediately close the pressure pot and apply 20 psi (1.4 bars) or greater pressure.
Note: Delay in applying pressure may introduce porosity.
4. Cure for 30 minutes.
5. Cool flask and remove denture from

the mold.

FINISHING AND POLISHING

Finish and polish the denture using normal procedures.

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York, PA 17401 U.S.A.
EC REP DeguDent GmbH